


P- 2- 22 Strip Sensor Pad 2" x 30" Mounts to Bed Edge / Rails etc.

USER MANUAL QUICK START GUIDE

Components:

1. Sensor Pad with Cord and RJ22 Connection.



Step 1.

Positioning the strip is the most important aspect of this unit. It could be placed along the edge of the bed to alarm when the patient swings their legs over to exit the bed. It should be placed from waist level down towards the feet. Place pad directly against the vinyl mattress using the Velcro tabs provided. Peel off the Velcro and mount opposite side from label. This unit provides no discomfort to the patient who is free to move around in the bed without alarming the unit. Uncoil the cord if not done already and run it straight down and underneath the mattress towards the back of the bed. One could also place the strip across the bed under the hips.

The strip might also be placed on a bed rail along the side of the bed to alarm when a patient begins to exit by grasping the rail. For a high risk patient the strip could be placed across the bed at a level with the patients back. When the patient begins to lift up to begin a process of getting out of bed the alarm will send a signal to the caregiver.



The RJ22 connection at the end of the cord may be connected into any of Emergency Caller Products alarm units such as the Patient Alert, the Extended Range Transmitter, TR-300, or the Long Range Transmitter, the LT-900. Consult your Emergency Caller Products sales representative to learn how these other products work with this model sensor pad. The pad may also be connected into other products that will receive the RJ22 connector.

Sensor Pads detect both Applied Weight "Patient ON" pad or Removed Weight "Patient Off". Selection made by switch on Universal Transmitter Unit or other monitoring products. Can be used to detect patient in bathtub or shower for entry or exit (weight off or weight on). Sensor pads are slip resistant and completely waterproof. Complete with velcro tabs and 8 ft. leads (except Chair-seat Pad with a 4 ft. lead). Lifetime sensor pads replace all disposable sensor pads. Emergency Caller Sensor pads have been used for over 25 years in security applications. They are durable, safe and withstand both outdoor and indoor environments. With lifetime guarantee you never again need to be concerned about constant replacement and expensive costs.