

Emergency Caller Products GETS HELP FAST!

(T) 800-227-2474 (F) 800-762-7613

CE-900 Wireless Enunciator QUICK START GUIDE


**Note* For use with
Emergency Caller Products Wireless
900MHz Supervised Transmitters Only!**



- 1 Components:
1. Enunciator (Receiver Horn-Strobe)
 2. Power Adapter
 3. Four (4) Velcro tabs



- 2 The Care Caller system is designed to assist the caregiver in monitoring patient activity. It is selectable for LOW or HIGH horn sound and comes with a Flashing Strobe. Selectable Off position when not in use. Unit has Aux. out (2.5mm jack, normally open) for Automatic Voice Dialer other connections.



- 3 Keep in mind the Enunciator has two sound levels- low and high. Also remember the transmitter range will be dependent to some degree on the structure of the residence. Also, you will need to place the unit within reach of an electrical outlet in order to plug in the adapter.



The Enunciator may be placed at Nurse Station, desk, on counter, nightstand, wall mounted or outside (covered). The caregiver should also consider where he or she might frequently be to be within sound of the enunciator when it alarms.

Important! Be sure to test your system often, with transmitter at locations where patient may be.



Now that we have decided where to place the enunciator let's go on to the next step. Begin by plugging the adapter into the top of the Enunciator by PWR as shown.


Next, plug the power adapter into the wall socket.



Now you are ready to test the Enunciator. Move the switch on the side of the unit to either low or high, for the sound output.

Test the system.

Activate the transmitter and listen for the sound. If you do not hear the Enunciator sound, check to be sure the switch is in Lo or Hi position and that the adapter is plugged in the wall. When your Enunciator sounds you have successfully setup and tested your system.



Reset Button

When activated the transmitter will send a signal to the Enunciator and the alarm sounds and the strobe flashes. The Enunciator **continues to sound until reset button is pressed** by Caretaker, assuring Help Has Arrived! Hold the button down for at least ten (10) seconds. The sound is then silenced, until transmitter is activated again. The Enunciator may be turned off when operation not desired.

Congratulations! Your system is all set up. Now you can feel Safe, Secure and at Ease whenever patients, family members, or friends are confined to bed, home or facility and need fast assistance. Care Caller sends the signal to call for help up to 1000 feet indoors (depending on structure, metal obstacles etc.). Be sure to test your system often as range may be effected by conditions.

Note: The enunciator emits beeping sound whenever Emergency Caller Transmitter battery is getting low.

See Reverse side for a list of other ECP Products.

ALL ECP wireless receivers have the capability to supervise all ECP universal transmitters. They also have the ability to pickup low battery warning signals from these transmitters. Every 60 seconds the transmitter sends an 'I'm here' signal to the receiver. If not sensed by the receiver in a thirty minute period the alarm will sound. This means protection from low power, non-working transmitters or transmitters removed from the facility.

Order Information: Ready-to-Go Complete Systems or Customize your own

CC-900-Care Caller-Enunciator with Pendant Transmitter. wrist band, necklace and belt clip.

Fall Alert Complete Systems - 900 MHZ Long Range - includes accessories

- FA-900- 1 Bed/Rail System with LT-900 transmitter, Enunciator, 1" x 30" Sensor Strip Pad, AC Adapter
- FA-900- 2 Bed/Rail System with LT-900 transmitter, Enunciator, 2" x 30" Sensor Strip Pad, AC Adapter
- FA-900- 3 Wheelchair-Seat System with LT-900 transmitter, Enunciator, 9"x16" Sensor Pad, AC Adapter
- FA-900- 4 Bed System with LT-900 transmitter, Enunciator, 14" x 24" Sensor Pad, AC Adapter
- FA-900- 5 Bed System with LT-900 transmitter, Enunciator, 18" x 24" Sensor Pad, AC Adapter
- FA-900- 6 Floor System with LT-900 transmitter, Enunciator, 24" x 30" Sensor Pad, AC Adapter
- FA-900- 7 Floor System with LT-900 transmitter, Enunciator, 24" x 36" Sensor Pad, AC Adapter

Emergency Caller Products presents advanced technology products to address those who need dependable care!



You will be very happy with the NEW ERU-900. Whenever there's an Emergency, a quick press of the Pendant (wear it on your wrist, neck, or hand-held) sends the signal to "Get Help Fast!" The Caller is triggered to deliver the emergency message! The Pendant is water resistant and can even be worn in the shower or bathtub! There are so many wonderful features with the

Wireless Fall Prevention Systems and Products including New Automatic "Man Down" transmitters that eliminate push buttons and pull cords. Also models with Universal transmitter that feature "pause", "sound alarm and transmit signal" when activated. Transmitters operate with Lifetime Sensor Pads and other sensors. Wireless Receiver Controllers include local PERS units, Enunciator and

WORKS AUTOMATICALLY! Detects patient wandering away and then sounds the alarm and calls for help! The Pendant Transmitter automatically sends a signal every 10 seconds to the Receiver-Dialer. If the receiver doesn't hear from the pendant within a 3 minute period, (indicating resident has traveled a distance away) the receiver-dialer will automatically sound the Alarm Strobe and/or dial emergency numbers and deliver the emergency message. No Monitoring Fees!

Lifetime Sensor Pads are available in many styles. Wheelchair (seat), bed, and floor. Exclusive Sensor Strip for beds and rails, allow for total comfort and freedom of movement in bed. Pads are thin and resident comfortable. Sensor pads are slip resistant and Completely Waterproof! Standard with velcro mounting tabs and 8 ft lead (chair-seat pads with 4 ft lead). Pads quick lock connect to transmitters and devices with RJ-22 plug. Lifetime Sensor Pads DIRECTLY REPLACE all disposable sensor pads. With Lifetime Guarantees, you'll never again need con-



Emergency Caller Products (T) 800-227-2474 (F) 800-762-7613
Made in USA