

**1.** The wireless pendant is attached to the patient with a non-removable band and sends a continuous signal to the Wanderer Alert unit every 10 seconds.



**2.** When the patient attempts to wander away, the Wanderer Alert unit detects the missing signals and sounds the horn/strobe alarm and/or dials the telephone, cell or pager numbers that you have installed and leaves your message in any language. A switch on unit allows you to select dial out.

**3.** The Wanderer Alert system is capable of monitoring both High Risk and Low Risk patients.

## WPT-900 Wireless Pendant